PREPARATORY STEPS

To ensure **consistent** testing conditions and **effective** tracking of body composition changes, follow these guidelines before each **InBody Test**.



Maintain normal fluid intake the day before and hydrate one hour prior to testing



Test at the same time of the day every time



Use the restroom prior to testing May skew fat mass/PBF result



Stand upright for 5-10 minutes prior to testing



Allow 20min for body temp to stabilize after being exposed to cold/hot temperature May skew fat mass/PBF result

DO NOT... $_{ extstyle -}$



Exercise prior to testing

May skew fat mass/PBF or SMM result Test only after fully recovered (1-3 days)



Eat prior to testing

May skew fat mass/PBF result Must wait at least 3 hours after last meal



Wear heavy clothing, accessories, or jewelry while testing May skew fat mass/PBF result



Consume alcohol or caffeine 24 hours prior to testing



Have lotion/ointment on hands or feet May skew fat free mass result



Wear any socks or pantyhose while testing May skew fat mass/PBF result

TESTING CONSIDERATIONS



Cosmetic Implants May skew fat mass/PBF result



Low carb/Keto Diet May skew fat mass/PBF result



Metal Implants May skew fat free mass result

Any skewed results will be consistent from test to test, so we recommend using your initial test as a baseline for tracking your body composition over time.



